



**FRESENIUS
KABI**

caring for life

KabiAccess

Non-dedicated Sets

The KabiAccess portfolio of non-dedicated sets includes extension sets, gravity sets, and add-on filters. It provides the most clinically convenient, safe and efficient configurations for your IV systems.



Short Extension Sets

FK20200E

6" (153mm) Minibore extension set with slide clamp, male luer lock with non-bonded needleless connector

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	6	0.3	100



FK20201E

6" (153mm) Macrobore extension set with pinch clamp male luer lock with non-bonded needleless connector

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	6	1.2	100



FK20203E

8" (20.32cm) Macrobore High Pressure extension set with purple pinch clamp, male luer lock with bonded needleless connector, 325 PSI


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	8	1.2	100



FK20210E

5" (130mm) Bi-fuse Minibore extension set with 3 slide clamps male luer lock and 2 non-bonded needleless connectors


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	5	0.5	100



FK20211E

5" (130mm) Bi-fuse Minibore extension set with 2 slide clamps male luer lock and 2 nonbonded needleless connectors

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	5	0.5	100




Short Extension Sets

FK20215E

5" (130mm) Tri-fuse Minibore extension set with 3 slide clamps male luer lock and 3 non-bonded needleless connectors

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	5	0.8	100



Extension Sets

FK20220E

18" (458mm) Macrobore extension set with roller clamp, 1 needleless connector, male luer-lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	2.4	100



FK20225E

18" (458mm) Macrobore extension set with roller clamp, male luer-lock and female luer


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	2.4	100



FK20245E

32" (813mm) Macrobore extension set with slide clamp, 2 needleless Y-site, male luer lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	32	5.0	48




Filter Extension Sets

FK20251E

18" (458mm) Microbore extension set with 0.22 Micron Filter pinch clamp, male lure-lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	4.6	100



FK20255E

18" (458mm) Macrobre extension set with 1.2 Micron Filter, pinch clamp, male lure-lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	4.6	100




Gravity Sets

FK30100G

100" (2540mm) macrobore standard gravity set with universal spike, 20 drops, roller clamp, slide clamp, male spin luer lock


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	100	18.0	48



FK30106G

106" (2692mm) Macrobre Primary Gravity set with universal spike, 20 drops, back check valve (BCV), 2 needleless Y-sites 90" and 10" above distal end (patient), slide clamp, roller clamp, male spin luer lock

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	106	19.0	48




Syringe Pump Extension Sets

FK50030S

30" (762mm) Minibore extension set with slide clamp, male luer-lock and female luer


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	30	0.5	100



FK50060S

60" (1524mm) Minibore extension set with slide clamp male luer-lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	60	1.0	100



For more information contact your sales representative or call customer service at 888-386-1300 or visit: www.agiliasystem.com

Copyright © 2021 Fresenius Kabi AG. All rights reserved.
All trademarks shown are the property of their respective owners.

Refer to the Instructions for use for a complete list of warnings and precautions associated with the use of this product.



**FRESENIUS
KABI**
caring for life

Fresenius Kabi USA
Three Corporate Drive
Lake Zurich, Illinois 60047
855-FK-INFUSION
855-354-6387
www.fresenius-kabi.us/agilia

0101-KA-01-05/21