

KabiAccess

Non-dedicated Sets

The KabiAccess portfolio of non-dedicated sets includes extension sets, gravity sets, transfusion sets, syringe pump extension sets, and secondary sets. It is your **ally** for convenient and efficient configurations for your IV systems.




Short Extension Sets

FK20200KE

6" Minibore extension set with slide clamp, male luer lock with non-bonded K-Zero® needleless connector


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.3	6	0.039	50



FK20201KE

6" Macrobore extension set with pinch clamp, male luer lock with non-bonded K-Zero needleless connector


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1.2	6	0.114	50



FK20202KE

6" Macrobore extension set with pinch clamp, male luer lock with non-bonded K-Zero needleless connector (Non PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1.2	6	0.114	50



FK20203KE

8" Minibore extension set with pinch clamp, male luer lock with K-Zero needleless connector, high pressure


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	8	0.039	50



FK20204KE

9.5" Minibore bi-fuse extension set pressure rated, 2 K-Zero needle-free connectors, 2 slide clamps, spin male luer lock, high pressure


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.4	9.5	0.039	50



FK20205KE

10" Macrobore bi-fuse extension set, 2 K-Zero needleless connectors, 2 slide clamps, spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2	10	0.114	50



FK20210KE

5" Bi-fuse Minibore extension set with 3 slide clamps, male luer lock, 2 non-bonded K-Zero needleless connectors

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	5	0.039	50




The diagram shows a 5-inch bi-fuse minibore extension set. It features a male luer lock connector on the left, followed by a section of tubing with three slide clamps. The tubing then splits into two parallel lines, each ending in a non-bonded K-Zero needleless connector.

FK20211KE

5" Bi-fuse Minibore extension set with 2 slide clamps, male luer lock, and 2 K-Zero needleless connectors

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	5	0.039	50




The diagram shows a 5-inch bi-fuse minibore extension set. It features a male luer lock connector on the left, followed by a section of tubing with two slide clamps. The tubing then splits into two parallel lines, each ending in a K-Zero needleless connector.

FK20215KE

5" Tri-fuse Minibore extension set with 3 slide clamps, male luer lock and 3 non-bonded K-Zero needleless connectors

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.7	5	0.039	50




The diagram shows a 5-inch tri-fuse minibore extension set. It features a male luer lock connector on the left, followed by a section of tubing with three slide clamps. The tubing then splits into three parallel lines, each ending in a non-bonded K-Zero needleless connector.

Extension Sets

FK20220KE

18" Macrobore extension set with roller clamp, 1 K-Zero needleless connector, male luer lock and female luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2.4	18	0.114	50




The diagram shows a 18-inch macrobore extension set. It features a male luer lock connector on the left, followed by a section of tubing with a roller clamp. The tubing then ends in a female luer lock connector and a K-Zero needleless connector.

FK20225KE

18" Macrobore extension set, roller clamp, male luer lock and female luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2.4	18	0.114	50




The diagram shows a 18-inch macrobore extension set. It features a male luer lock connector on the left, followed by a section of tubing with a roller clamp. The tubing then ends in a female luer lock connector.

FK20226KE

20" Macrobore extension set with 2 K-Zero needleless Y-sites and male luer lock adapter with retractable collar

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
3.3	20	0.114	50




The diagram shows a 20-inch macrobore extension set. It features a male luer lock adapter with a retractable collar on the left, followed by a section of tubing with two K-Zero needleless Y-sites. The tubing then ends in a female luer lock connector.

Long Extension Sets

FK20245KE

35" Macrobre extension set with slide clamp, 2 K-Zero needleless Y-sites, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
5	35	0.114	50



FK20246KE

40" Extension set with 3-way macrobre stopcock with rotating male luer lock, 2 K-Zero needleless Y-sites, and male luer lock with retractable collar

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
6.3	40	0.114	50




Filter Extension Sets

FK20251KE

18" Macrobre filter extension set with 0.2 micron filter, pinch clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



FK20252KE

18" Macrobre filter extension set with 0.2 micron filter, pinch clamp, male luer lock and female luer lock (Non-PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



FK20255KE

18" Macrobre filter extension set with 1.2 micron filter, pinch clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



FK20256KE

18" Macrobre filter extension set with 1.2 micron filter, pinch clamp, male luer lock and female luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50




Gravity Sets

FK30100KG

100" Universal 20 drop spike, roller clamp, slide clamp, male spin luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
---------------------	-----------------------	------------------	------------

18 100 0.114 50




FK30104KG

104" Universal 10 drop spike, back check valve (BCV), roller clamp, 3 K-Zero needleless Y-sites (4", 23" and 83") from 2-piece male luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
---------------------	-----------------------	------------------	------------

15 104 0.114 50




FK30106KG

106" Universal 20 drop spike, back check valve (BCV), 2 K-Zero needleless Y-sites (10" and 90") above distal end, slide clamp, roller clamp, and male spin luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
---------------------	-----------------------	------------------	------------

19 106 0.114 50




FK30110KG

110" Universal 60 drop spike, back check valve (BCV), 2 K-Zero needleless Y-sites (10" and 90") above distal end, slide clamp, roller clamp, and male spin luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
---------------------	-----------------------	------------------	------------

19 100 0.114 50



Transfusion Sets

FK40080KT

80" 2 bag spikes, 2 roller clamps, 20 drop drip chamber with 200 micron filter, roller clamp, slide clamp, 1 K-Zero needleless Y-Site 10" above distal end, male spin luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
---------------------	-----------------------	------------------	------------

25 80 0.114 50




FK40100KT

100" 2 bag spikes, 2 roller clamps, 20 drop drip chamber with 200 micron filter, hand pump, roller clamp, slide clamp, 1 K-Zero needleless Y-Site 10" above distal end, male spin luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
---------------------	-----------------------	------------------	------------

45 100 0.118 25




Syringe Pump Extension Sets

FK50030KS

30" Minibore extension set with slide clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	30	0.039	50



FK50032KS

30" Macrobores extension set with slide clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2	30	0.114	50



FK50034KS

30" Microbore extension set with slide clamp, male luer lock and female luer lock (Non-PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	30	0.039	50



FK50060KS

60" Minibore extension set with slide clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1	60	0.039	50



FK50062KS

60" Macrobores extension set with slide clamps, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
3	60	0.114	50



FK50064KS

60" Microbore extension set with slide clamp, male luer lock and female luer lock (Non-PVC Fluid Path)

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1	60	0.039	50






Secondary Sets

FK60034KSC

34" Macrobre secondary set and extension hook

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
10	34	0.114	50



FK60035KSC

34" Macrobre secondary set

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
10	34	0.114	50



Please refer to the instructions for use for a complete list of safety and caution information for KabiAccess Non-Dedicated Administration Sets.

KabiAccess Non-Dedicated Sets are not made with natural rubber latex or DEHP.
All measurements are approximate.

Your ally in transforming infusion delivery



IVENIX
INFUSION SYSTEM

Engineered from the ground up to leverage the power of the digital age.



agilia
Connect

Agilia Connect integrated system of pumps, software, sets, accessories and services.



freeflex[®]
IV Bags & Solutions

A growing portfolio of IV solutions in **freeflex** and **freeflex+** IV bags.



Fresenius Kabi USA **Nutrition**

Our parenteral nutrition products are used to help care for critically ill and chronically ill patients in hospitals, long-term care facilities and at home.



For more information, contact your IV Therapy Sales Representative at **1.855.354.6387** or visit www.kabiaccess.com.

All trademarks shown are the property of their respective owners.